#### PRECISE CONDUCTIVITY METER





#### HIGHLIGHTS:

- Modern and functional housing
- Excellent price-performance ratio
- o 3-line display/overhead display at the push of a button
- Backlit
- Waterproof (IP67)
- o Robust, long battery life
- o including high-quality measuring cells for a wide range of applications
- o quick determination of measured values



#### G 1410

Item No. 610006

Universal conductivity meter up to max 100 mS/cm, including graphite measuring cell

#### G 1410-1002

Item No. 474039

Universal conductivity meter Device, measuring cell LF 202, 2-pole graphite, permanently connected, in case GKK1002

#### G 1420

Item No. 610007

high-resolution ultrapure water conductivity meter up to max 100  $\mu$ S/cm, including stainless steel measuring cell

### G 1420-1002

Item No. 474040

Display:

high-resolution ultrapure water conductivity meter Device, measuring cell LF 200 RW, 2-pole stainless steel/PEEK, firmly connected, in case GKK 1002

#### GENERAL:

When developing the new G 1000 series, the focus was consciously on the essential functions of the measurement technology. Pure measurement with a focus on precision, speed and reliability, packaged in a compact housing, enables the outstanding price/performance ratio, Made in Germany

Made in Germany.
The new handheld meters also feature an ergonomic design, IP 65/67 dust and water protection, and an illuminated display. The compact conductivity meter, the G 1410, has a precise and long-lasting wide-range measuring cell for universal use from DI water to seawater. As the G 1420, it has a specialized measuring cell for high-resolution pure/ultrapure water applications.

# APPLICATION:

Fresh and sea water aquariums, reverse osmosis and similar filters, cleaning processes, cooling/lubricating processes, plant cultivation and agriculture; Laboratories, quality assurance, service

	quality distributes, set vice	
TECHNICAL SPECIFICATIONS:	G 1410	G 1420
	Wide range measuring instrument including graphite measuring cell	Ultrapure water version including stainless steel measuring cell
Measurement:	Conductivity, salinity, TDS	Conductivity, specific resistance
Measuring ranges:	each with automatic measuring range switching	
Conductivity:	02000 µS/cm 0.0020.00 mS/cm 0.0100.0 mS/cm	0.0002.000 μS/cm 0.0020.00 μS/cm 0.0100.0 μS/cm
Specific resistance:	-	0.01000.2000 MOhm * cm 0.0102.000 MOhm * cm 0.0120.00 MOhm * cm
TDS:	02000 mg/l	
Salinity (PSU):	0.050.0 g/kg	
Temperature:	-5.0 + 105.0 °C	-5.0 + 105.0 °C
accuracy		
Conductivity:	$\pm$ 0.5 % of meas. MW. $\pm$ 0.5 % FS	Type. $\pm$ 1 % of meas. MW. $\pm$ 0.5 % FS
Temperature:	± 0.3 °C	± 0.3 °C
Temperature compensation:	off: deactivated nLF: non-linear, according to EN 27888	off: deactivated nLF: non-linear, according to EN 27888 LIN: linear with adjustable coefficient NaCl: For weak NaCl solutions according to EN 60746-3
Reference temperatures:	20 and 25 ℃	20 and 25 ℃
Sensors/measuring inputs:	permanently connected 2-pole measuring cell with integrated temperature sensor	
Measuring cell:	2-pole measuring cell, Ø 12 mm (graphite), cable 1.2 m (others at extra cost)	2-pole measuring cell, Ø 12 mm (stainless steel 1.4404, 1.4435), cable 1.2 m (others at extra cost)

Service: 4 durable, easy-to-use buttons Additional automatic measuring range switching, functions: automatic temperature compensation Working Device: -20 .. + 50 °C, 0..95 % RH (non-condensing) conditions: Measuring cell: -5 .. + 80 °C (briefly 100 °C) 2 x AA batteries, > 1000 h battery life Power supply: Protection class: IP65/IP67 Housing: break-proof ABS housing 108 x 54 x 28 mm (H x W x D) without Dimensions: sensor connection approx. 200 g (G 1410) Weight: approx. 230 g (G 1420) Device with measuring cell, test report, Scope of delivery: 2 x batteries, operating instructions

# ACCESSORIES OR SPARE PARTS:

#### **GKL 100**

Item No. 601396

Conductivity control solution Control solution 1413  $\mu$ s/cm, 100 ml bottle,

according to DIN EN 27888

### **GKL 101**

Item No. 601398

Conductivity control solution 84  $\mu s$  /cm, 250 ml bottle

# GKL 102

Item No. 601400

Conductivity control solution 50 ms/cm, 100 ml bottle

#### HD-22-3

Item No. 700040

Freely positionable laboratory sensor holding arm for sensors Ø12mm

### GWZ-01

Item No. 603499

Flow-through chamber for sensors with Ø 12 mm, tube connection Ø 6 mm

### ST-G1000

Item No. 611373 Protective bag, leather

#### GB AA

Item No: 610049

Replacement battery Mignon (AA) 1.5 V (2 pieces are required)

### **GKK 1002**

Item No. 411907

Case for G1000 series water analysis small

3-line with battery level indicator, backlit, protected by break-resistant glass, overhead display at the push of a button